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
None

Reporting Code Changes

The following manufacturers have either deleted or updated codes for this mailing.

None

Storage and Stability Instructions

1. Specimens are stable for 7 days when stored unopened at 2 - 8°C.
2. Once open, specimens must be tested within 5 days.

Critical Reporting Information

1. The following portions of this Survey are formally evaluated: Anti-A Titer Result (Exact), and Anti-D Titer Result (Exact).
2. For any testing that you do not routinely perform in your laboratory, leave all reporting areas for that test blank unless otherwise noted.
3. The method used in performing this Survey should be the one routinely used in your laboratory. While laboratories are encouraged to consider using the uniform procedure, results using the uniform procedure should be reported only if it is the procedure you routinely use.
4. Please record on the result form whether you used the "uniform procedure" or "other procedure."
5. Information about the uniform procedure can be found in the ABT-A 2008 Part I Summary Report and for reference from the result form, see article: Aufbucher, Widi-Egger J, Dunst L, for the Biomedical Excellence for Safer Transfusion Initiative; the Transfusion Medicine Resource Committee of the College of American Pathologists. Reducing the variation in performance of antibody titration. *Transfus Lab Med* 2008;132:1194-1201.

Important: See the Backstage Warning at the end of these instructions.
Customer Contact Center 800-323-4040 option 1 (domestic), or 847-832-7900 option 1 (international)

ABT-A 2013 (ART-ABT)

KITS 25403889 9 01 01

CAP# Y136472-01 SEQ# 01 CD 04 0123

PROJECTS 001ART 1 ART 2

LABOR 100
WEST CHINA MEDICAL UNIVERSITY HOSP
751400 28-0000000

LABOR 28-0000000

Page: 1
Results must be received at the CAP no later than
midnight, Central Time by the due date below:
July 02, 2013

Online: cap.org (preferred method)
Fax: 866-FAUCAP (866-365-0027)

Antibody Titer Survey Result Form

Reporting Code Selection

- If your method summary page states, "Please Provide a
"KID CODE,"
 - If your code is listed incorrectly, or
 - If you have changed your methodology.
- Review the master list for all appropriate codes and enter it on the result form. If there is no master list, select the code directly on the result form.
- If you cannot find an appropriate code
- Select Other from the kit instructions or result form and describe your method in the line of Other section of the result form.
- If you need assistance, please call the Customer Contact Center at 800-323-4040 option 1 (domestic), or 847-832-7900 option 1 (international).
- View the e-LAB Solutions™ user guide via cap.org

Anti-A Titer -- ABT-01 with titer cell ABT-03

Summary of Uniform Procedures*

Dilution of Sample

Volumes: 1 mL
Diluent: 0.9% NaCl Solution (Normal Saline)
Technique: Coloured paper, new for each dilution

Tube Technique: Anti-A

First reading:
Incubation: Room temperature for 30 minutes and read without additional agglutination.

Second reading

Incubation: 15°C for 30 minutes
Washes: Four times with at least ten-fold the volume of red cell + sample
Testing phase: Anti-IgG
Endpoints: no, read macroscopically

Gel Card Technique: Anti-A

Two separate cards
Gel cards: incubated at room temperature for 15 minutes
Anti-IgG: incubated at 37°C for 15 minutes
Endpoint: no

Provision of Reagent Red Cell

For the purposes of ABT Survey testing, participants should use the titer cell provided with the ABT kit.

*Note: Any step not specified above is left to the definition of the laboratory.

Customer Contact Center 800-323-4040 option 1 (domestic),
or 847-832-7900 option 1 (international)


54085

APN

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J-A 2014

Kit Instructions
PCAP 2014

Comprehensive Transfusion Medicine, Educational Challenge, and Electronic Crossmatch Survey

J • JE1 • EXM

Table of Contents

Kit Contents	1
Important: Before You Begin	1
Detailed Testing Instructions	2
Reporting Your Results	4
Behavioral Warning	5
For Assistance	5
Master Link	7
ABC Group/Rh Type Worksheet	8

Kit Contents

J	J-OR1 — J-OR17 J-OR18 — J-OR19
JE1	JE-OR1 and JE-OR2
EXM	EXM-OR1 — EXM-OR3 simulated donor units
	EXM-OR4 — EXM-OR8 corresponding red cells

Important: Before You Begin

Now for this mailing

Supplemental questions have been added to the result form. The Transfusion Medicine Resource Committee encourages all participants to volunteer responses to these questions. Results will be reported in the Participant Summary report.

Reporting Code Changes

The following manufacturers have either deleted or updated codes for this mailing. None.


Storage and Stability Instructions

- Store vials at 2 - 8°C.
- The red blood cell suspensions should be tested within 14 days of receipt. Perform testing as soon as possible.

Critical Reporting Information


- Antibody Identification:** If your laboratory performs Antibody Identification, you must complete the Antibody Identification section of the result form, even if the screen was negative.
- Compatibility Testing:**
 - Perform and record serologic results even if an ABO/Rh incompatibility is expected. Incompatible crossmatches are included in the Survey to verify that your laboratory is proficient in detecting them.
 - For the purpose of this proficiency testing, serologic crossmatches must be performed.

Important: See the Behavioral Warning at the end of these instructions.
Customer Contact Center 800-323-4040 option 1 (domestic) or 801-847-832-7000 option 1 (international)

J-A
2014

OR

KITS 27189437 3 01 30

CAP# T156X72-01 BRQ# 01 01 64 0348

PROBTS: 2014 J A 1

SHANGHAI
WEST CHINA SECOND UNIVERSITY HOSP
TEL#86-20-85021541 FAX#86-20-85555888

Page: 1

Results must be received at the CAP no later than
mailing. Central Time by the due date below.

March 12, 2014

Online: cap.org (preferred method)
Fax: 800-FAX-2CAP (800-323-2227)

Comprehensive Transfusion Medicine Survey Result Form

Reporting Code Selection

- If your method summary page states, "Please Provide a Valid Code."
 - If your code is listed incorrectly, or
 - If you have changed your methodology.
- Review the master list for an appropriate code and enter it on the result form. If there is no master list, select the code directly on the result form.
- If you cannot find an appropriate code:
Select Other from the list instructions or result form and describe your method in the Use of Other section of the result form.
If you need assistance, please call the Customer Contact Center at 800-323-4040 option 1 (domestic) or 801-847-832-7000 option 1 (international).

View the e-LAB Solutions™ user guide via cap.org

© CAP 2014


ABO Group

Manufacturer Code	Anti-A	Anti-B	Anti-A,B	A,B Cells	B Cells
010	010	010	010	010	010
29	29	26	26	29	29
Exception Code	☐ 11 ☐ 22				
ABO Group J-018/J-015	ABO Group J-020/J-025	ABO Group J-030/J-035	ABO Group J-040/J-045	ABO Group J-050/J-055	
☒ 100 Group A ☐ 101 Group A+ ☐ 102 Group A- ☐ 103 Group B ☐ 104 Group B+ ☐ 105 Group B- ☐ 106 Group AB ☐ 107 Group AB+ ☐ 108 Group AB- ☐ 109 Group O ☐ 110 Screen/Phenotype and cell group do not agree; additional testing or sample required	☐ 111 Group A ☐ 112 Group A+ ☐ 113 Group A- ☐ 114 Group B ☐ 115 Group B+ ☐ 116 Group B- ☐ 117 Group AB ☐ 118 Group AB+ ☐ 119 Group AB- ☐ 120 Group O ☐ 121 Screen/Phenotype and cell group do not agree; additional testing or sample required	☐ 122 Group A ☐ 123 Group A+ ☐ 124 Group A- ☐ 125 Group B ☐ 126 Group B+ ☐ 127 Group B- ☐ 128 Group AB ☐ 129 Group AB+ ☐ 130 Group AB- ☐ 131 Group O ☐ 132 Screen/Phenotype and cell group do not agree; additional testing or sample required	☐ 133 Group A ☐ 134 Group A+ ☐ 135 Group A- ☐ 136 Group B ☐ 137 Group B+ ☐ 138 Group B- ☐ 139 Group AB ☐ 140 Group AB+ ☐ 141 Group AB- ☐ 142 Group O ☐ 143 Screen/Phenotype and cell group do not agree; additional testing or sample required	☐ 144 Group A ☐ 145 Group A+ ☐ 146 Group A- ☐ 147 Group B ☐ 148 Group B+ ☐ 149 Group B- ☐ 150 Group AB ☐ 151 Group AB+ ☐ 152 Group AB- ☐ 153 Group O ☐ 154 Screen/Phenotype and cell group do not agree; additional testing or sample required	
Ungraded	Ungraded	Ungraded	Ungraded	Ungraded	Ungraded
155 Anti-A identified	156 Anti-A identified	157 Anti-A identified	158 Anti-A identified	159 Anti-A identified	

Customer Contact Center 800-323-4040 option 1 (domestic)
or 801-847-832-7000 option 1 (international)

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